



Be Right™



Flask, Filtering, Polypropylene, 1000 mL

Kat. číslo produktu:	2089953
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Polypropylene filtering flask with side-arm at 45° angle prevents tipping. Barbed side-arm accepts 1/4-3/8" I.D. tubing for application of vacuum. See page 380 for tubing; see page 375 for stopper.

Technické údaje

Jednotky:	1 each
Kapacita:	1000 mL
Materiál skrine:	Polypropylene