



Be Right™



Free chlorine indicator solution for chlorine analyser CL17/ CL17sc (473 mL)

Kat. číslo produktu:

2314011

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dostupné



Technické údaje

| | |
|------------------------------------|---|
| Doba použiteľnosti: | 5 months |
| Parameter: | Chlorine, Free |
| Platforma: | Roztok |
| Počet testov: | 450 |
| Podmienky skladovania: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |
| Prístroj: | CL17/CL17sc |
| Rozsah merania: | N/A |
| Volume / package size e Post Text: | 473 mL |