



PHG301 General-purpose pH electrode, glass

Kat. číslo produktu:	E11M001
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

pH glass electrode with Ag/AgCl reference system. Measuring range 0-12 pH, temperature range 0-80°C. Min. immersion depth 12mm.

Technické údaje

Dĺžka:	103 mm
Materiál:	Glass
Obsah balenia:	Includes: Radiometer Analytical PHG301-9 pH electrode with Certificate of Conformity.
Rozsah merania:	0 - 12 pH pH
Teplotný rozsah:	0 - 80 °C
Teplotný senzor:	none
Typ elektrody:	Laboratory pH electrode (sensor only, external reference probe needed)
Typ elektródy:	PHG301
Typ senzora:	Glass
Záruka:	6 mesiacov

Obsah balenia

Includes: Radiometer Analytical PHG301-9 pH electrode with Certificate of Conformity.

Potrebné príslušenstvo

- CL114 cable, 1m, for screw cap electrode FX/S7/coax, BNC plug (Item A94L114)
- REF321 Reference electrode, Ag/AgCl (Item E21M002)