




Batoh na prenosné prístroje, veľký, nylon, s kolieskami

Kat. číslo produktu:	2985000
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Pre pohodlnú prepravu prenosných prístrojov a až 6 štandardných roztokov 50 ml.

Praktický

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Prispôsobiteľný

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

Ergonomický

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

Technické údaje

Dĺžka:	152 mm
Hmotnosť:	Kapacita: 13.6 kg
Materiál:	Nylon, umelé vlákno
Rozmery (V x Š x H):	432 mm x 356 mm x 152 mm
Výška:	ID: