



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Kat. číslo produktu:	2842100
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Technické údaje

Materiál:	Glass Microfiber
Priemer:	47 mm mm
Veľkosť balenia:	100
Veľkosť póru:	1 µm