



## ISE25Cl-9 Chloride ion selective electrode (ISE) half-cell, screw cap

**Kat. číslo produktu:**

**E41M003**

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

### **The right probe for your application.**

To meet your analysis needs, the Radiometer Analytical ion selective electrode range offers high quality sensors for a variety of ionic species in solution.

### **Factory controlled to ensure high quality and premium performance.**

Each Radiometer Analytical Ion Selective Electrode is carefully controlled in production. Each ISE is delivered with an individual control report including information on the response curve, checkpoint value, and slope.

---

## **Technické údaje**

Záruka:

6 mesiacov