



Addista Multiparametrový štandardný roztok

Kat. číslo produktu:	LCA700
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Systém na zabezpečenie analytickej kvality (AQA), ktorý sa skladá z jedného kombinovaného štandardného/spajkovacieho roztoku a dvoch roztokov na medzilaboratórne porovnanie skúšok.

Použiteľné pre kyvetové testy

- LCK304 amoniak, 0,015 - 2,0 mg/L $\text{NH}_4\text{-N}$
- LCK311 chlorid, 1 - 70 mg/L Cl^-
- LCK228 draslík, 5 - 50 mg/L K
- LCK328 draslík, 8 - 50 mg/L K
- LCK348 ortofosforečnany, 0,5 - 5,0 mg/L $\text{PO}_4\text{-P}$
- LCK414 / LCK1414 CHSK, 5 - 60 mg/L O_2
- LCK238 Laton, 5 - 40 mg/L TN_b

Confidence in your results

Continuous documentation of Analytical Quality Assurance measures using the Addista system substantiates the correctness of your water testing results.

Less complexity

Each Addista standard solution contains several parameters. This reduces the number of bottles with single parameter standards in your lab significantly.

Less error sources

The Addista solutions are ready-to use with Hach's LCK cuvette test. Since there is no dilution needed this allows less hands-on-time and eliminates an additional error source when analysing standard solutions.

Technické údaje

Platforma:	LCA
Podmienky skladovania:	2 - 8 °C (uchovávať v chlade)
Popis:	Multiparametrový štandardný roztok ADDISTA