



Vrečkový kolorimeter II Súprava kolorimetra pre stanovenie chloridov kyvetovým testom

Kat. číslo produktu: 5953000V.07

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Portable Colorimeter programmed for the determination of Chloride with cuvette tests. In carrying case, with batteries and operating instructions.

Following Hach cuvette test can be evaluated with this instrument:

LCK311, 1-70 mg/L Cl / 70-1000 mg/L Cl

Cuvette test has to be ordered separately.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to lightning and large figures on the display

Reliable results with no mains connection

Rugged construction yet lightweight

Technické údaje

Parameter:	Chloridy (vhodné pre LCK311)
Rozsah merania:	1 - 70 mg/L / 70-1000 mg/L Cl