



Beaker, polypropylene, low form, 50 mL

Kat. číslo produktu:	108041
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technické údaje

Jednotky:	1 each
Kapacita:	50 mL
Materiál:	Polypropylene, Low Form
Materiál skrine:	Polypropylene
Priemer, vonkajší:	51.1 mm with spout
Veľkosť balenia:	1 each