



Beaker, polypropylene, low form, 600 mL

Kat. číslo produktu:	108052
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Technické údaje

Číslo metódy:	8006
Jednotky:	1 each
Kapacita:	600 mL
Materiál:	Polypropylene, Low Form
Priemer, vonkajší:	110.7 mm with spout
Prístroj:	All except DR/4000, PC II