



Be Right™



Funnel, analytical, 254 mL approx. Volume

Kat. číslo produktu:

108370

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

Technické údaje

Veľkosť balenia:

1 each