



Buret, Acrylic, 50 mL

Kat. číslo produktu:	1405941
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Durable, crystal-clear acrylic body is individually calibrated to meet the accuracy requirements of ASTM E287. Leak-proof stopcock with self-lubricating PTFE TFE plug. Excellent for general lab or field use in schools, industrial, municipal or food applications. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Can contain liquids within a pH range of 1 - 14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents.

Technické údaje

Kapacita:	50 mL
Materiál skrine:	Acrylic
Špeciálne funkcie:	Buret
Stupnice:	0.1 mL