



Test kit, nitrogen, nitrate color disc, model NI-14

Kat. číslo produktu:

1416100

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model:	NI-14
Model testovacej súpravy:	NI-14**
Názov metódy:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate
Platforma:	Farebný disk
Počet testov:	50
Rozsah merania:	Dual Range
Rozsah merania (2):	0 - 10 mg/L NO ₃ ⁻ -N
Smallest increment:	0.02 0.2
Štandardný rozsah merania:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Štýl kufríka:	B