



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Kat. číslo produktu:	190902
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Technické údaje

Body Material:	Borosilicate Glass
Kapacita:	60 mL
Materiál:	Borosilicate Glass
Package size UOM:	1pk
Veľkosť balenia:	1 each