



Test kit tannin/lignin

Kat. číslo produktu:

193701

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

Tannin and lignin are analyzed together and reported as mg/L tannic acid.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	Tannin/Lignin
Model testovacej súpravy:	TA-3
Názov metódy:	Color disc/Tyrosine
Parameter:	Tannin/Lignin
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0 - 15 mg/L as Tannic Acid
Rozsah merania (2):	0 - 150 mg/L
Smallest increment:	0.5
	5
Štandardný rozsah merania:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
Štýl kufríka:	R