



Cone, Imhoff settling, 1 L, polycarbonate, w/ screw closure

Kat. číslo produktu:	206700
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Fully graduated cones for determination of small amounts of settleable matter in sewage by Imhoff Sedimentation Method. Glass cone features reinforced tip for ruggedness, while the Polycarbonate cone features a leak-proof, removable closure for withdrawal of sediment for analysis. Support ring, support stand, and cone holders available and sold separately.

Technické údaje

Kapacita:	1 L
Veľkosť balenia:	1 each