



Test kit, chloride, model CD-51

Kat. číslo produktu:

208601

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	CD-51
Názov metódy:	Drop count titration/ Silver nitrate
Parameter:	Chloride - As Cl^-
Platforma:	Drop Count Titration
Počet testov:	100
Poznámka:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	500 - 10000 mg/L & 5000 - 100000 mg/L Cl^-
Smallest increment:	500 5000
Štýl kufríka:	B