



Be Right™



Bottle, Sampling, Low Density, Polyethylene, 1 L, 6/pk

Kat. číslo produktu:	2087083
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Technické údaje

Body Material:	Low Density Plastic
Kapacita:	1000 mL
Materiál:	LDPE
Package size UOM:	6pk
Veľkosť balenia:	6 /pk