



Cylinder, Graduated, Mixing, Glass, 10 mL ± 0.1 mL, 0.2 mL divisions, glass stopper #13

Kat. číslo produktu: 2088638
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Technické údaje

Certifikácia:	Class B
Kapacita:	10 mL
Materiál skrine:	Glass
Metóda merania:	To Contain
Typ:	S jednou stupnicou, zmiešavací
Veľkosť balenia:	1 each