



Flask, Erlenmeyer, polymethylpentene capacity 125 mL

Kat. číslo produktu:	2089843
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Technické údaje

Jednotky:	1 each
Kapacita:	125 mL
Materiál:	Polymethylpentene
Typ:	Erlenmeyer