



Crucible, Gooch, 25 mL Capacity

Kat. číslo produktu:	2098640
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Perforated bottom, glazed inside and out except for outside bottom surface. Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Technické údaje

Kapacita:	25 mL
Veľkosť balenia:	1 /pk