



Crucible, Gooch, 40 mL Capacity

Kat. číslo produktu:	2098641
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Technické údaje

Kapacita:	40 mL
Veľkosť balenia:	1 /pk