



Test kit, formaldehyde, model FM-1

Kat. číslo produktu:

2183100

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	FM-1
Názov metódy:	Drop count titration/ Thymolphthalein
Parameter:	Formaldehyde - As CH_2O
Platforma:	Drop Count Titration
Počet testov:	100
Rozsah merania:	0,05 - 1% & 0,5 - 10%
Smallest increment:	0.05% 0.5%
Štýl kufríka:	D