



Bag, sterile, WhirlPak, 207 mL, 500/pk

Kat. číslo produktu:	2233100
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Approved by U.S. Environmental Protection Agency (USEPA). Disposable, transparent polyethylene bags. Gas-sterilized, pull-tab openings. Can be marked with waterproof felt pen or marker. Sure "whirl" closure. Nasco.

Technické údaje

Kapacita:	207 mL
Veľkosť balenia:	500