



Test kit, nitrogen , model NI-8

Kat. číslo produktu:

224100

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model testovacej súpravy:	NI-8
Názov metódy:	Color disc/Nessler reagent
Parameter:	Nitrogen, Ammonia, mid range - For freshwater
Platforma:	Farebný disk
Počet testov:	100
Prepravná hmotnosť:	1
Rozsah merania:	0-3,0 mg/L NH ₃ -N
Smallest increment:	0.1
Štýl kufrika:	B