



Be Right™



Agar R2A, dehydrated 500 g

Kat. číslo produktu:	2281034
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Low nutrient medium for enumeration and cultivation of bacteria from potable water. Used for pour or spread plate and MF methods; described in the 16th edition of Standard Methods. Recommended for predicting nitrification more accurately in potable systems using chloramines as a secondary disinfectant.

Technické údaje

Media Type:	R2A Agar
Veľkosť balenia:	500 g