



Chlorid sodný standardní roztok, 18 000 $\mu\text{S}/\text{cm}$, 500 mL

Kat. číslo produktu: 2307449
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

Conductivity standard solutions for laboratory and field analyses

Standard konduktivity. 18 000 $\mu\text{S}/\text{cm}$ (10 246 mg/L NaCl). 500 mL.

Technické údaje

Koncentrácia:	18000 $\mu\text{S}/\text{cm}$ at 25°C $\pm$ 50 $\mu\text{S}/\text{cm}$ at 25°C
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Roztok
Podmienky skladovania:	10 - 25 °C (protect from light)
Popis:	Conductivity calibration solution
Veľkosť balenia:	500 mL
Vzorec:	NaCl