



Be Right™



Test kit, ascorbic acid, model ASC-1

Kat. číslo produkta:

2308100

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:

ASC-1

Názov metódy:

Drop count titration/

Iodometric

Parameter:

Ascorbic Acid

Platforma:

Drop Count Titration

Počet testov:

100

Poznámka:

*mg/L unless otherwise noted, ppb = $\mu\text{g}/\text{L}$, ppm = mg/L; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Rozsah merania:

10 - 200 mg/L

Smallest increment:

10

Štýl kufríka:

D