



A-1 broth medium, dehydrated, 500 g

Kat. číslo produktu:

2309934

Nebezpečné

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.



A-1 Broth Medium, Dehydrated. Broth for 24-hour MPN method for fecal coliforms. See Microbiological Testing for prepared medium. Bottle/500 g.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Media Type: Živné médium A-1

Region: Globální

Veľkosť balenia: 500 g