



Be Right™



Butterfield's buffered phosphate

Kat. číslo produktu:

2319110

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

Butterfields Buffered Phosphate, 450 mL bottle. pk/10.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Jednotky:	10 ks v balení
Objem / Veľkosť balenia:	450 mL
Veľkosť balenia:	10/pk