



Filter, acrocap, pore size 0.2 µm

Kat. číslo produktu:	2323810
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Technické údaje

Filter Type:	Filtr, AkroCap
Packing:	NA
Priemer:	NA
Veľkosť póru:	0,2 µm