



Test kit, molybdate, model MO-LR

Kat. číslo produktu:

2359300

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Technické údaje

Model:	MO-LR
Model testovacej súpravy:	MO-LR
Názov metódy:	Color disc/ Ternary complex
Parameter:	Molybdate
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	Single Range
Smallest increment:	0.25
Štandardný rozsah merania:	0 - 3.0 mg/L Mo
Štýl kufríka:	D