



Filter, membrane, pore size 0.2 µm, diameter 25 mm, 100/pk

Kat. číslo produktu: 2387200

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Filters are treated with an irgalan black process that renders the filters dark gray, eliminates natural autofluorescence and aids visibility. These polycarbonate screen membrane filters exhibit extremely low non-specific adsorption. The fluorescing dye used does not discolor the filter. Manufacturing provides exceptionally clean membranes, as free as possible of viable and nonviable microorganisms. Diameter: 25 mm. Pore size: 0.2 µm. Poretics.

Technické údaje

Filter Type:	Filtr
Materiál:	Polycarbonate
Packing:	100 ks v balení
Priemer:	25 mm
Veľkosť balenia:	100 /pk
Veľkosť póru:	0,2 µm