




Filter, acrodisc, pore size 0.2 µm

Kat. číslo produktu:	2415600
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Technické údaje

Filter Type:	Filter, acrodisc
Materiál:	Polysulfone
Packing:	NA
Priemer:	NA
Veľkosť póru:	0.2 µm