



## Ampules, m-TGE/TTC broth, plastic, pk/50

Kat. číslo produktu:	2428450
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

m-TGE w/TTC, 2 mL plastic ampules. pk/50.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Technické údaje

Materiál:	Plastic
Media Type:	m-TGE/TTC Broth
Veľkosť balenia:	50 50/pk