



Be Right™



Flask, Erlenmeyer, Wide Mouth, 2000 mL, Subdivisions 600 mL

Kat. číslo produktu:	2489454
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included.

Technické údaje

Jednotky:	1 each
Kapacita:	2000 mL
Materiál:	Glass Wide Mouth