



Stirring rods, hand stirrer, disposable nylon, 10/pk

Kat. číslo produktu:	2496510
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Hand-held mini-mixer used with test tubes, micro centrifuge tubes and containers too narrow for conventional mixers. Stirring rod flares out for use in larger spaces. Clear vinyl carrying case, 11.4 x 6 cm (4.5 x .2.4") stirring rod, 2 AA batteries included. Dimensions including rod 20.3 x 3.2 x 2.5 cm (8 x 1.25 x 1"). Weight 3 oz.

Technické údaje

Veľkosť balenia: 10 /pk