



Sample container, 45 mL

Kat. číslo produktu:	2552114
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Plastic, leak-proof. Corning Snap Seal™ offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Technické údaje

Kapacita:	45 mL
Priemer, vonkajší:	30 mm (approx)
Veľkosť balenia:	600 /pk