



## Sample container, 120 mL, pk/250

|                      |                 |
|----------------------|-----------------|
| Kat. číslo produktu: | 2552212         |
| EUR Cena bez DPH:    | Kontaktujte nás |
| Dodávka do 2 týždňov |                 |

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

---

### Technické údaje

|                    |                 |
|--------------------|-----------------|
| Kapacita:          | 120 mL - 250/pk |
| Priemer, vonkajší: | 45 mm (approx)  |
| Veľkosť balenia:   | 250 /pk         |