



Test kit, surface water

Kat. číslo produktu:

2559800

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal Tester and all necessary apparatus and reagents for testing.

Ideal tool for both industrial and municipal water treatment plants

Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature

Rugged, lightweight kit includes the convenient pH Pocket Pal[™] Tester and all necessary apparatus and reagents

Technické údaje

Aplikácie:	Surface Water
Platforma:	Farebný disk
Počet testov:	100