



Agar, tryptic soy, dehydrated 100 g

Kat. číslo produktu: 2565926

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

General purpose medium used without blood for cultivating fastidious microorganisms.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Technické údaje

Media Type:	Tryptic Soy Agar
Veľkosť balenia:	100 g