



Filter, TSS glass fiber, pore size 1.5 µm, diameter 110 mm, 100/pk

Kat. číslo produktu:	2569100
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Technické údaje

Materiál:	Glass Microfiber
Priemer:	110 mm
Veľkosť balenia:	100 /pk
Veľkosť póru:	1.5 µm