



**Be Right™**



## Filter assembly, portable, complete (millipore microfil)

Kat. číslo produkta: 2586200

### ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Navrhované náhrady:

- 2586201

For membrane filtration features unique membrane lifting device. Depressing the lever breaks the vacuum and lifts the membrane above the 47 mm membrane support for easy retrieval with sterilized forceps. Designed for use with Microfil funnels (see Filtration Funnels), creating an efficient, low-maintenance system for membrane filtration in the lab or field. Rugged metal construction with vacuum side-arm at the base. Works with any vacuum source.

---

### Technické údaje

Veľkosť balenia:

1