



Be Right™



Scoop, sampler, sterile

Kat. číslo produktu:

2600900

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Sterile scoop samplers hold various capacities. One-piece molded scoops used to obtain samples with one scoop. White, high-impact styrene. Individually bagged, lengths measuring from 17.8 to 21.6 cm (7 to 8.5"); 100 to case. Capacity: 60 mL (20 oz). Bowl: 25 x 76 x 38 mm (1 x 3 x 1.5") Bel-Art.

Technické údaje

Kapacita: 60 mL

Veľkosť balenia: 1