



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Kat. číslo produktu:	2613049
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Technické údaje

Farba:	Transparent
Jednotka:	6 pk
Kapacita:	500 mL
Materiál:	LDPE
Stem:	Blue Cap/Stem
Veľkosť balenia:	6 pk