



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Kat. číslo produktu:	2618800
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Technické údaje

Filter Type:	Filter
Materiál:	Polyethersulfone
Priemer:	47 mm
Veľkosť balenia:	200 /pk
Veľkosť póru:	0.2 µm