



Sample vial, 120 mL w/out thiosulfate, 200/pk

Kat. číslo produktu:	2634002
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Technické údaje

Dĺžka:	93 mm
Kapacita:	Fill line at 100 mL ($\pm$ 2.5 mL)
Priemer:	43 mm
Špeciálne funkcie:	Sample Vial, without Thiosulfate
Veľkosť balenia:	200/pk