



Cylinder, Graduated, 1000 mL, Certified

Kat. číslo produktu:	2636253
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Technické údaje

Certifikácia:	Class A
Kapacita:	1000 mL
Materiál skrine:	Glass
Metóda merania:	To Deliver
Typ:	S jednou stupnicou
Veľkosť balenia:	1 each
Výška:	465 mm