



Cylinder, graduated, mixing, certified, 25 mL

Kat. číslo produktu:	2636340
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Technické údaje

Certifikácia:	Class A
Kapacita:	25 mL
Materiál skrine:	Glass
Metóda merania:	To Contain
Typ:	S jednou stupnicou, zmiešavací
Veľkosť balenia:	1 each