



Filter Unit, Type A

Kat. číslo produktu:

2656600

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Certified for water quality work, this sterile, disposable, polystyrene filter unit contains a 47 mm, 0.45 µm cellulose nitrate membrane for superior recovery and growth of microorganisms. It offers good chemical resistance to dilute acids, aromatic and aliphatic hydrocarbons, and conforms to American Public Health Association (APHA) and Standard Methods. Filter unit holder is available for use in multiple filtrations.

Technické údaje

Priemer: 47 mm

Veľkosť balenia: 12 per pack